

GASTRO WORKSTATION

GI Diagnostic Function Testing System



Alpine's Gastro Workstation is the complete solution for all your GI diagnostic testing needs.

The workstation provides flexibility by allowing you to expand to multiple applications while offering the consistency and reliability of working with the same software platform.

Applications:

- Ambulatory pH testing
- Esophageal manometry
- Anorectal manometry
- Biofeedback
- Sphincter of Oddi manometry
- Small bowel manometry
- Electrogastrography (EGG)

THE GASTRO WORKSTATION IS AVAILABLE IN THREE VERSIONS

Regional versions available with applicable language and power supply.

STANDING VERSION

- High cart 110 cm (43.5") high
- Large storage drawer
- Large locking anti-static wheels
- Dimensions: 22.5" x 22" x 43.5" (57.2 cm x 55.9 cm x 110.5 cm)
- Includes proprietary CPU, R/W CD drive, 19" flat screen color monitor, built-in hospital-grade transformer, and color inkjet printer. (US only)
- Alpine Biomed PC unit

SITTING VERSION

- Low cart 80 cm (31.5") high
- Large locking anti-static wheels
- Dimensions: 21.6" x 22.0" x 31.5" (55 cm x 56 cm x 80 cm)
- Includes proprietary CPU, R/W CD drive, 19" flat screen color monitor, built-in hospital-grade transformer, and color inkjet printer. (US only)
- Alpine Biomed PC unit

NOTEBOOK VERSION

- Portable clamshell design
- Includes R/W CD drive and color inkjet printer (US only)

SOFTWARE

- GastroTrac® software disk
- Microsoft® Windows® XP Professional, SP2
- Microsoft® Internet Information Services
- Microsoft® Internet Explorer 6.0
- Optional remote access for Technical Support (US only)

Alpine Biomed Corp.

17800 Newhope St., Suite B

Fountain Valley, CA 92708

USA

Tel: +1 800 222 0074

Tel: +1 714 839 8426

Fax: +1 714 839 8429

www.alpinebiomed.com

Alpine Biomed ApS

Tonsbakken 16-18


DK - 2740 Skovlunde

Denmark

Tel: +45 44 57 90 00

Fax: +45 44 57 90 10

www.alpinebiomed.com



REQUIRES

- Digitrapper® pH400 for pH testing
- Polygraf™ ID multi-parametric physiologic recorder for Manometry, EGG and Biofeedback
- GastroTrac® licenses to perform desired tests
- Accessories

ADDITIONAL HARDWARE / ACCESSORIES AVAILABLE

- Digitrapper® pH Monitoring System
- IR dongle for pH400 upload
- GastroTrac® licenses for advanced pH, esophageal, anorectal, Sphincter of Oddi, small bowel Manometry, EGG, and Biofeedback
- Water-perfused, solid-state (US only) and air-charged manometry catheters, as well as other manometry accessories
- Manometric Perfusion Pump (US only)
- Shelf for Manometric Perfusion Pump
- Biofeedback and EGG accessories
- pH catheters
- Calibration buffer for pH catheters
- Inkjet printer
- Catheter Puller

Alpine Biomed and logo, Digitrapper, and GastroTrac are registered trademarks and/or registered service marks of Alpine Biomed or one of its subsidiaries, in the United States and/or other countries. Polygraf is a trademark of Alpine Biomed. The absence of a name or logo in this list does not constitute a waiver of any intellectual property rights that Alpine Biomed or its subsidiaries have established in any of their products, features, or service names or logos. Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries.

Disclosure: See product labeling for specific directions regarding product use, contraindications, warnings, precautions and potential adverse events.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

F061207B
©Alpine Biomed Corp. 2009
All Rights Reserved